

THANK YOU...

Thank you for investing in your independence with HealthCraft Products. We are confident that you will find the unique design and durable construction of this product to be of great assistance in your daily routine. To fully enjoy your e2™ Ceiling Mount Trapeze, please take a few moments to read these instructions. You should store these instructions in a safe place for future reference.

SAFETY CONSIDERATIONS

CAUTION: It is your responsibility to see that your e2™ Ceiling Mount Trapeze is properly assembled, installed, and cared for. Failure to follow instructions in this document could result in serious injury.

e2™ CEILING MOUNT TRAPEZE INSTALLATION / MAINTENANCE

The e2™ Ceiling Mount Trapeze should only be used in locations where the ceiling surfaces are structurally sound. Recommended weight limit is 300 lbs. No more than 150 lbs should be exerted on trapeze handle (sitting up assistance only). If you are not equipped to undertake the outlined work, we would recommend that you have your e2™ Ceiling Mount Trapeze installed by a qualified contractor. Inspect, or have inspected (and tighten if necessary) the screws on a monthly basis.

WARRANTY

All HealthCraft™ Products are covered by a one year limited warranty with the exception of hand grips which are subject to normal wear. Buyer hereby indemnifies, agrees to hold harmless and defend HealthCraft™ Products Inc. from and against any and all liabilities, claims, (founded and unfounded), losses, damages, costs and expenses (including without limitation consequential damages and reasonable professional fees) resulting from buyers specification, application, or improper use of goods described hereon; buyers omission or neglect.

HealthCraft™ Products Inc. does not assume any liability for damage resulting from services performed by others or faulty installation, misuse or misapplication of goods sold by HealthCraft™ Products Inc.

HealthCraft™ Products Inc. shall not be liable for prospective profits or special, indirect, or consequential damages, or for the cost of any corrective work done without HealthCraft™ Products Inc. prior approval. HealthCraft™ Products Inc. total liability hereunder shall in no event exceed the purchase price of the goods specified hereon.

INSTALLATION INSTRUCTIONS

1. Determine the location of the grab bar as shown in **Figure 1**. The actual user of the e2™ trapeze should lie down in bed in their normal resting position with their hands placed palms down on their thighs. The grab bar will be mounted such that the strap hangs directly down over the palms of the user as shown by "X". Hold the trapeze handle under where the grab bar will go, with the upper rail of the handle closest to the head of the bed. Ensure that the trapeze handle is not too close to allow a comfortable sit up action. The grab bar should also be installed so that it extends from the centre of where the user will lie down to an area that they can slide the trapeze when not in use, as shown in "Z" in **Figure 1**.

2. Once the location is determined, mark the hole placement of the grab bar on ceiling. It is necessary that both grab bar flanges are secured to joists in the ceiling as shown in **Figure 2**.

3. Adjust the length of the strap to set the appropriate height of the trapeze (when ones fingers can comfortably grasp around the bottom rail of the trapeze when lying down in bed as shown by "Y" in **Figure 1**). This is done by pulling some slack through the cross pieces of the strap buckle and pulling the buckle down to shorten the strap, or sliding it up to lengthen the strap.

4. Open the loop on the trapeze strap and slide trapeze handle on.

5. Place the hook on the ceiling grab bar with the open side facing the foot of the bed, as shown in **Figure 1**.

PARTS DIAGRAM




Figure 1.




Figure 2.

