


Assista-Rail™ with Assista-Cane™

A standard rail for home beds with accessories to meet changing needs. The Assista-Cane™ is a fixed handle that attaches to the Assista-Rail™ to improve standing and transfer support.

Features and Benefits

- **Easy Installation** - Rail installs in minutes with no tools. Simply slide the rail between mattress and box spring. Includes safety strap and anti-slide grips to hold rail secure.
- **Sleep Safely** - Effective rail design allows for in-bed mobility & bed queuing (i.e. gives an idea of where the mattress edge is located).
- **Transfer Safely** - By minimizing upper body twisting and reducing wrist strain, the optional Assista-Cane™ provides additional support to reduce the risk of a forward fall.
- **Dual Rail System** - Install rails on both sides of the bed for added support (conversion kit available).



Technical Specifications:

Assista-Rail™	Heavy gauge steel, 7/8" / 2.25 cm diameter. High temperature baked powder coat gloss white enamel. Extends under mattress 24" / 61 cm, up from the mattress 16" / 41 cm and along the side of the bed 18" / 46 cm.
Grip Material	Grey, NPVC non-latex grip. 27" / 68.5 cm long x 1.25" / 3 cm diameter.
Traction Grip	Grey, premium quality, non-latex grip. 4.5" / 11 cm long x .84" / 2 cm diameter.
Safety Strap	Strap 1" / 2.5 cm x 100" / 254 cm length with side release buckle.
Assista-Cane™	Heavy gauge steel, 7/8" / 2.25 cm diameter. High temperature baked powder coat gloss white enamel. 31" / 79 cm tall x 7.75" / 20 cm, with grey PVC end cap & rubber cane tip.
Grip Material	Grey, NPVC non-latex grip. 11.5" / 29 cm long x 1.25" / 3 cm diameter.
Weight Capacity	250lbs / 113 kg.
Assistance Level	



ENTRAPMENT WARNING: Due to the potential risks associated with head / limb entrapment, the Assista-Rail™ is not recommended for users with the following conditions: erratic body movements, altered mental status, pre-existing conditions such as confusion, restlessness and dementia.