



# Advantage Pole™ / Portable Options

A floor mounted pole. Excellent where ceiling restrictions prevent a floor to ceiling pole. Perfect beside a toilet, bath, bed, or chair.

## Features and Benefits

### • Safe, Effective Support.

Push or pull on the pole to assist with sitting, standing and moving.

### Reduce Caregiver Strain.

- Maximize the use of available strength while a caregiver supervises. Promotes safe and independent movement.

### • Portable Model (sold separately).

Pole quickly secures and releases from a low profile, floor mounted plate. Install additional floor plates to use pole in multiple locations.



## Technical Specifications:

Advantage Pole™	54" / 137 cm x 1.5" / 3.81 cm diameter. White powder coated steel.
Grip Length	20" / 51 cm. Grey latex free grip.
Advantage Pole™ Mast	28" / 71 cm height x 2.25" / 5.7 cm diameter.
Floor Plate	7.5" / 19cm diameter floor plate (screw pattern 6 x 3/8" holes).
Weight Capacity	Bariatric weight capacity 450 lbs / 205 kg.