



Assista-Cane™

A quarter length bed assist rail that allows for a fixed cane to be added as needs change. Provides a useful combination of bed safety with bed mobility. Fits twin to king size beds.

Features and Benefits

▪ **Simple Installation.**

Slide the Assista-Rail™ between the box spring and mattress, secure to bed frame with mandatory safety strap. Traction grips help to prevent rail from shifting during use.

▪ **Assista-Cane™ Add On Kit (sold separately).**

A sturdy handle that mounts to the Assista-Rail™ and can be fixed at any angle to provide a safe weight bearing surface. Minimizes the risk of falls by reducing body twisting.

▪ **Assista-Cane™ Leg Extender (sold separately).**

Increases the Assista-Cane™ height from 32.5" / 83 cm to 35.5" / 90 cm for taller users.

▪ **Dual Rail System (sold separately).**

Conversion kit allows an Assista-Rail™ to be placed on both sides of a bed.



Assista-Cane™ Specifications:

Assista-Cane™	(LxW): 32.5" x 8" / 83 cm x 20 cm. 7/8" / 2.22 cm diameter. Powder coated (white) steel.
Grip	11.625" / 29.5 cm x 1.25" / 3 cm diameter. Grey, latex free grip.
Weight Capacity	300lbs / 136 kg.

Assista-Rail™ Specifications:

Assista-Rail™	(LxWxH): 33.5" x 20" x 17.25" / 85 cm x 51 cm x 44 cm. 7/8" / 2.22 cm diameter. White powder coated steel.
Grip	27" / 68.5 cm x 1.25" / 3 cm diameter. Grey, latex free grip.
Traction Grips	4.5" / 11 cm. Grey, latex free grip (x4).
Safety Strap	1" / 2.5 cm x 100" / 254 cm length strap with side release buckle.
Weight Capacity	300lbs / 136 kg.