



SuperPole™ Horizontal Rail System

A compact horizontal rail system based on the innovative SuperPole™. Central to the system are floor to ceiling poles, that provide a secure mounting structure for the horizontal rails.

Features and Benefits

- **Walking Therapy.** Horizontal Rail System can help to improve recovery time, regain mobility and independence. Excellent guide for short walking distances. Available as a Parallel Bar or Single Rail System.
- **Install And Relocate In Minutes.** Easy setup and removal. Ideal for situations where a permanent installation may not be an option. Great for home or professional use.
- **Height Adjustable Rails.** Secure at virtually any height. An extremely low range, 18" / 46 cm from floor, makes it ideal for pediatric applications. Loosen the captive knobs to remove rails in seconds.



Parallel Bar System



Single Rail System

Technical Specifications:

| | |
|------------------------|---|
| Horizontal Bar | 100" / 254 cm x 1.5" / 3.81 cm diameter. White powder coated steel. |
| Height Adjustable | 18-40" / 46-102 cm from floor. |
| SuperPole™ | 90" / 229 cm x 1.5" / 3.81 cm diameter. White powder coated steel. |
| Grip Length | 20" / 51 cm. Washable, latex free grip. |
| Floor to Ceiling Range | 93-99" / 236-253 cm (standard). |
| Expansion Mechanism | 1-1/8" fine cut thread. Zinc plated steel jackscrew. |
| Top Plate | 20" / 51 cm length. Top plate rubber will not mark ceiling. |
| Base Plate | 5" / 13 cm diameter. Non-marking, base plate rubber ribbed for better gripping. |
| Weight Capacity | 300 lbs / 136 kg. Bariatric weight capacity 550 lbs / 249 kg. |

SuperPole™ Accessories (each sold separately):

| | |
|-------------------------|---|
| Uni-Fit Extender™ | Increase floor to ceiling range from 100-120" / 254-305 cm. Not to be used with the SuperTrapeze™, Bariatric, & Ultra Models. |
| SuperPole™ Ultra | Longer pole accommodates floor to ceiling range from 121-140" / 307-356 cm. Not to be used with the SuperBar™, Safety Island™ or SuperTrapeze™. |
| Ceiling Plate Extender™ | Allows standard ceiling plate 20" / 51 cm length to span 30" / 76 cm with holes for ceiling attachment to 24" / 61 cm ceiling beam centers. |