



ModPole™

A revolutionary bed mounted support rail. For moderate support getting into and out of bed. Perfect for home or hospital beds.

Features and Benefits

- **Safe.** Push or pull on the pole to assist with sitting, standing, and safe repositioning. Will not restrict getting into and out of bed.
- **Reduce Caregiver Strain.** Maximize the use of available strength while a caregiver supervises. Promotes safe and independent movement.
- **Easy Installation.** Mounts securely to angle iron bed frames of box spring and mattress home beds. Also available for Invacare® and Hillrom® hospital beds.

NOTE: For home beds, mounting bracket includes a telescopic beam that spans 23-40" / 58-102 cm. Fits most single to king size bed frames. See website for details.



Technical Specifications:

ModPole™	36" / 91 cm x 1.5" / 3.8 cm diameter.
Grip Length	20" / 51 cm. Grey latex free grip.
Height Adjustable	4" / 10 cm, overall height range dependent on bed height.
Weight Capacity	325 lbs / 147 kg.