



ModRail™

A revolutionary bed mounted support rail. A fixed dual bed rail (ModRail™) with a pivoting horizontal bar. For moderate support getting into and out of bed. Perfect for home or hospital beds.

Features and Benefits

- **Simple To Use.** A gentle lift unlocks the bar. Keep the bar up to move freely. Let the bar down to lock. Provides 8 locking positions. Bar is height adjustable.
- **Safe.** Push or pull on the bar to assist with sitting, standing, and safe repositioning. To help prevent a forward fall, the bar will lock if released. Bed rail helps define bed edge and can help reduce falls while sleeping. Will not restrict getting into and out of bed.
- **Reduce Caregiver Strain.** Maximize the use of available strength while a caregiver supervises. Promotes safe and independent movement.
- **Easy Installation.** Mounts securely to angle iron bed frames of box spring and mattress home beds. Also available for Invacare® and Hillrom® hospital beds.

NOTE: For home beds, mounting bracket includes a telescopic beam that spans 23-40" / 58-102 cm. Fits most single to king size bed frames. See website for details.



Technical Specifications:

ModRail™	14" / 36 cm x 1.25" / 3.2 cm diameter. Locks every 45° (8 positions). 4" / 10 cm spacing between rails on side gate. Latex free grip.
Top Handle	12.5" / 32 cm (removable).
Height Adjustable	4" / 10 cm, overall height range dependent on bed height.
Weight Capacity	325 lbs / 147 kg.