



Stainless Steel SuperPole™

A professional grade, expansion fit, floor to ceiling pole. Features stainless steel construction designed to withstand frequent and rigorous use. A great choice for home or care facilities.

Features and Benefits

• Well Placed Support.

Install virtually anywhere support is needed, unlike grab bars which are limited to the nearest wall.

• Safe Bathing.

Ideal for accessible showers where well-placed support can help to reduce slips and falls. A 3-way seal in the base prevents standing water build up and the growth of bacteria.

• Simple Installation.

Hold the ceiling plate and turn the pole. Pole will expand between floor and ceiling. Lock in place with the provided set screws.

• Will Never Rust.

Stainless steel will never rust, stain or corrode. Provides years of worry free maintenance.

Intended for use with "built-in place" tile/fiberglass showers with 100% floor contact support (i.e. set in cement). Not recommended for bathtubs/shower enclosures raised from floor with gap spacing between the enclosure and the floor or ceiling.



Technical Specifications:

SuperPole™	90" / 229 cm x 1.5" / 3.81 cm diameter. Stainless steel.
Knurled Section	20" / 51 cm.
Floor to Ceiling Range	93-98" / 236-249 cm (standard).
Expansion Mechanism	.5" fine cut thread. Located at the top of the pole.
Top Plate	20" / 51 cm length. Top plate rubber will not mark ceiling.
Base Plate	5" / 13 cm diameter. Non-marking, base plate rubber ribbed for better gripping.
Weight Capacity	450 lbs / 205 kg.