

CAUTION
RISK OF INJURY!
 Understand these instructions prior to using product (Only for use on INVACARE #5411 Hospital Style Bed)
VERTICAL LOAD BEARING ONLY - Product intended for moderate VERTICAL weight bearing on unlocked rail only. - DO NOT use if sideward (lateral) support is needed.
AVOID RAIL LATCH UNLOCKING - Ensure that rail is at least 1"/25mm off floor when bed is in **lowest** position (the face plate moves down with the bed, and rail could become unlatched)
 Maximum user weight 250lbs/113kg
 Do not use on wet or slippery floors
 Product intended for moderate assist with sitting up and standing only.

