



Grab Bars




SuperPole


Easy Installation.

1-2hrs

locating studs, levelling,
drilling, sealing, etc.



15mins



It is recommended that grab bars are secured directly to structural wood studs with a minimum of 2 screw per flange. This requires locating studs and drilling holes through ceramic tile. If installed in the bathroom, it is important to ensure a watertight seal to avoid long term damage (from mold and mildew).

The SuperPole installs in minutes by simply turning a jackscrew at the bottom of the pole. This action raises the pole until the top plate securely contacts the ceiling - no drilling required. The SuperPole must be located directly under a ceiling joist for maximum support.

No Damages on Install.





Grab bars are secure, but they are permanent - not ideal for needs that may evolve (such as reselling a home, injury rehabilitation or changing health conditions). If the grab bars need to be removed, significant effort is required to repair the damaged tiles and holes.



The SuperPole installation is strong enough to support the user, yet it can be removed in minutes with no damage to the floor or ceiling. The top and bottom plates have rubber pads that provide traction and protect the surfaces from scratching. Remove or relocate the pole with ease.


Consistent Support.



There are many safety products for different areas of the home. Grab bars for the bathroom, bed rails for the bed, and some areas such as the sofa has few available products. These products all work differently and have various rail diameter, which can lead to confusion.

Using a consistent product for all areas of the house will help increase confidence. The SuperPole can be installed beside a bed, inside the bathroom or beside a favour couch - anywhere it is needed.


User Optimized Support.



Many safety products are attached to fixed surfaces, which places support at the installation point and not ideally placed for the user. (For example, a grab bar on the wall, or a bedrail to the bed).

The SuperPole can be installed virtually anywhere there is a floor and ceiling. It can be placed exactly where required, especially when safety needs to be in the middle of a room.

Evolution of Support.



Grab bars can provide a basic level of support. Unfortunately, they are unable to adapt as the user's support requirements change.

The SuperPole can be used to add on support accessories to help as user's needs evolve. These support accessories attach directly onto the pole to help the user maintain independence.

Grab Bars



SuperPole

- | | | | |
|---|---|---|---|
| ✗ | Installation requires drilling thru tiles and drywall | ✓ | Expansion-fit pole easily installs with no drilling |
| ✗ | Removal requires additional work to fix holes in wall | ✓ | Can be removed / relocated with no damage to floor or ceiling |
| ✗ | Different products are required for each room in the home | ✓ | SuperPoles can be used in any room with floor and ceiling |
| ✗ | Support location determined by product installation convenience | ✓ | Support location determined by where user needs it |
| ✗ | Cannot upgrade safety as needs evolve | ✓ | Support accessories can be installed as needs evolve |